

Bencox[®] Bipolar Cup



Bencox[®] Bipolar Cup



Optimal Fit


- A large array of component allows precise patient matching for optimal patient fit.
- Outer diameters from 38 to 60mm allow the surgeon to meet the unique challenges posed in revision and extended situations.
- The bipolar cup offers use with 22, 28mm femoral heads.

* Bipolar Cup Size / Head Size

38 / 22	39 / 22	40 / 22	41 / 22	42 / 22	43 / 28	44 / 28	45 / 28	46 / 28	47 / 28	48 / 28	49 / 28
											
50 / 28	50 / 28	52 / 28	53 / 28	54 / 28	55 / 28	56 / 28	57 / 28	58 / 28	59 / 28	60 / 28	
											

Wide Range of Motion

- Enhanced 68° ROM resists the impingement to the UHMWPE locking ring.



Proportional Eccentricity

- Multipolarity is achieved via multiple offsets within the system as cup sizes increase, the amount of eccentricity increases.
- Provides eccentricity that is proportional to the patient's anatomy.



Simple Locking Ring Mechanism


- Ring mechanism allows for easy assembly of the femoral head.



Enhanced Stability

- Superior pull-out strength resists the dislocation of the femoral head.

Testing Laboratory: ENDOLAB (No. 265.100430.10.1482)
Date: 26 May 2010



*Reference: Kummer, FJ, Jaffe WL, Lesaka K, Perez F, Bull Hosp Jt Dis. 2005;63(1-2):5-8

Ordering Information

Part No.	Description	ASTM F799 (forged)	
		Part No.	Description
01.03.038	Bipolar Cup 38	H2.A51.0038	Bipolar Cup 38
01.03.039	Bipolar Cup 39	H2.A51.0039	Bipolar Cup 39
01.03.040	Bipolar Cup 40	H2.A51.0040	Bipolar Cup 40
01.03.041	Bipolar Cup 41	H2.A51.0041	Bipolar Cup 41
01.03.042	Bipolar Cup 42	H2.A51.0042	Bipolar Cup 42
01.03.043	Bipolar Cup 43	H2.A51.0043	Bipolar Cup 43
01.03.044	Bipolar Cup 44	H2.A51.0044	Bipolar Cup 44
01.03.045	Bipolar Cup 45	H2.A51.0045	Bipolar Cup 45
01.03.046	Bipolar Cup 46	H2.A51.0046	Bipolar Cup 46
01.03.047	Bipolar Cup 47	H2.A51.0047	Bipolar Cup 47
01.03.048	Bipolar Cup 48	H2.A51.0048	Bipolar Cup 48
01.03.049	Bipolar Cup 49	H2.A51.0049	Bipolar Cup 49
01.03.050	Bipolar Cup 50	H2.A51.0050	Bipolar Cup 50
01.03.051	Bipolar Cup 51	H2.A51.0051	Bipolar Cup 51
01.03.052	Bipolar Cup 52	H2.A51.0052	Bipolar Cup 52
01.03.053	Bipolar Cup 53	H2.A51.0053	Bipolar Cup 53
01.03.054	Bipolar Cup 54	H2.A51.0054	Bipolar Cup 54
01.03.055	Bipolar Cup 55	H2.A51.0055	Bipolar Cup 55
01.03.056	Bipolar Cup 56*	H2.A51.0056	Bipolar Cup 56*
01.03.057	Bipolar Cup 57*	H2.A51.0057	Bipolar Cup 57*
01.03.058	Bipolar Cup 58*	H2.A51.0058	Bipolar Cup 58*
01.03.059	Bipolar Cup 59*	H2.A51.0059	Bipolar Cup 59*
01.03.060	Bipolar Cup 60*	H2.A51.0060	Bipolar Cup 60*

* Bipolar Cup 56, 57, 58, 59, 60 Custom Order (available upon request)

Manufacturing Plant 12, Yeongsanhong 1-gil, Ipjang-myeon, Sebuk-gu, Cheonan-si, Chungcheongnam-do, Korea 31056
Tel : 82-41-585-7114 Fax : 82-41-585-7113

Seoul office 8F, 11F, Chungho Tower, 483, Gangnam-daero, Seocho-gu, Seoul, Korea 06541 Tel : 82-2-3445-5492~5 Fax : 82-2-3445-5467,5497

Corentec America 60 Washington Street, Suite 202, Morristown, NJ 07960, USA Tel : 1-862-777-8320 Fax : 1-862-777-8319

EC Representative : Emergo Europe Molenstraat 15, 2513 BH, The Hague, The Netherlands Tel : 31-70-345-8570 Fax : 31-70-346-7299